



These are the decorations on the child size cuff, which is shown smaller than actual size. Decorations for the infant size cuff are similar but not identical.

The actual color of both the child and infant cuffs is medium blue.

World Standard for blood pressure

The Decorated Calibrated® V-Lok® Cuff

Our Decorated cuffs are made of the same quality materials as our standard cuff. Familiar figures of animals and toys have been added in the hope that they may diminish a child's anxiety and confusion during part of a physical exam.

As with all **Baumanometer® Calibrated® V-Lok® Cuffs** our Decorated Cuffs feature:

- The Calibrated® V-Lok® system which immediately and clearly indicates whether or not the correct size is being used.
- Genuine Velcro® hook and loop fasteners which are precisely positioned and evenly matched in holding strength to provide the optimal number of open-close cycles.
- Cuffs that clearly indicate the center of the bladder on the inner and outer surfaces for ease of application on either right or left limb.
- The cuff material is a tightly woven, urethane coated, heat set, Dacron™ polyester fabric.
- The design, materials and high quality workmanship of products made in the USA.

Cat No.	Description
1871D	Child Decorated Calibrated® V-Lok® Cuff
1872D	Infant Decorated Calibrated® V-Lok® Cuff
1881D	Child Complete Cuff and Bladder
1882D	Infant Complete Cuff and Bladder
1821D	Child Complete Inflation System
1822D	Infant Complete Inflation System

Latex-free Decorated items can be ordered by adding NL to the appropriate catalog number, i.e. 1881DNL.

W.A. Baum Co., Inc.

620 Oak Street, Copiague, NY 11726 USA
 Tel: 631-226-3940 Fax: 631-226-3969
 www.wabaum.com info@wabaum.com

**Since 1916 Originator and Maker of Accurate
Blood Pressure Apparatus Exclusively**

Baumanometer®