



DETAILS.

Attention to detail has earned W.A. Baum Co., Inc. its reputation for designing and producing World Standard Blood Pressure Equipment for more than eighty-eight years.

The details found in every Calibrated® V-Lok® Cuff assure correct fit, ease of application, accurate, repeatable blood pressure readings and long service life.

Accurate: The Calibrated® V-Lok® system clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb.

Durable: Baum cuffs are made of tightly woven, heat set Dacron™ polyester fabric. They are double stitched at critical points for long wear. Hook and loop fasteners are evenly matched in holding strength to provide the optimal number of open-close cycles.

Comfortable: Cuffs conform to the shape of the limb and have no hard, stiff edges.

Safe: Each cuff is treated with an effective anti-microbial agent and can be paired with non-latex components for use with patients who are sensitive to natural rubber latex.

Guaranteed: Latex and Non-Latex inflation bags, bulbs and tubing are guaranteed for five years.

The Calibrated® V-Lok® Cuff - a World Standard ... in every detail.

Available through better-known medical equipment dealers worldwide.

For more information on the complete line of Baumanometer® Calibrated® V-Lok® Cuffs, call toll-free 1-888-281-6061 or e-mail: info@wabaum.com

W.A. Baum Co. Inc.

SINCE 1916

ORIGINATOR AND MAKER OF
BLOOD PRESSURE APPARATUS EXCLUSIVELY

W.A. Baum Co., Inc., Copiague, NY 11726 USA

Ordering Information

Calibrated® V-Lok® Cuff

Aneroid and Decorated cuffs also available. Catalog numbers are listed below.

Cat. No.	Size	Aneroid	Decorated	Limb Circumference
● 1870	Adult	1870AC	--	25 cm - 35 cm
1871	Child	1871AC	1871D	18 cm - 26 cm
1872	Infant	1872AC	1872D	10 cm - 19 cm
1873	Newborn	1873AC	--	6 cm - 11 cm
1874	Thigh	1874AC	--	46 cm - 66 cm
1868	Large Arm	1868AC	--	33 cm - 47 cm

Calibrated® V-Lok® Cuffs & Inflation Bags

Includes inflation bag with two tubes. For single tube cuff and bag sets, see catalog pages 12-13.

Cat. No.	Non-Latex	Size	Aneroid	Non-Latex	Deco.	Non-latex
● 1880	1880NL	Adult	1880AC	1880ACNL	--	--
1881	1881NL	Child	1881AC	1881ACNL	1881D	1881DNL
1882	1882NL	Infant	1882AC	1882ACNL	1882D	1882DNL
1883	1883NL	Newborn	1883AC	1883ACNL	--	--
1884	1884NL	Thigh	1884AC	1884ACNL	--	--
1869	1869NL	Large Arm	1869AC	1869ACNL	--	--



Note: All Latex items are natural rubber Latex.

Calibrated® V-Lok® Complete Inflation Systems



Comprises Cuff, Latex or Non-Latex Inflation Bag, Latex or Non-Latex Bulb and Air-Flo® Control Valve. For Aneroid and Decorated Complete Inflation Systems, use catalog numbers below. NIBP Inflation Systems are also available, refer to pages 14-16 of our equipment catalog.

Cat. No.	Non-Latex	Size	Aneroid	Non-Latex	Deco.	Non-latex
● 1820	1820NL	Adult	1820AC	1820ACNL	--	--
1821	1821NL	Child	1821AC	1821ACNL	1821D	1821DNL
1822	1822NL	Infant	1822AC	1822ACNL	1822D	1822DNL
1823	1823NL	Newborn	1823AC	1823ACNL	--	--
1824	1824NL	Thigh	1824AC	1824ACNL	--	--
1825	1825NL	Large Arm	1825AC	1825ACNL	--	--

Note: All Latex items are natural rubber Latex.



Color Coded Cuff Sizes



Color Code	Cuff	Cuff & Bag	Infl. Sys.
G-Gray	1873G	1883G	1823G
T-Tan	1872T	1882T	1822T
DG-Dark Green	1871DG	1881DG	1821DG
MB-Medium Blue	1870MB	1880MB	1820MB
S-Scarlet	1868S	1869S	1825S
DB-Dark Brown	1874DB	1884DB	1824DB

W.A. Baum Co., Inc.
620 Oak Street
Copiague, NY 11726 USA

Tel: 631-226-3940
Fax: 631-226-3969
Toll Free: 1-888-281-6061

Web: www.wabaum.com
Email: info@wabaum.com

W.A. Baum Co. Inc.

SINCE 1916

ORIGINATOR AND MAKER OF BLOOD PRESSURE APPARATUS EXCLUSIVELY