

Blood Pressure Replacement Parts



1820

Calibrated® V-Lok® Inflation System


Cat. No.	Non Latex(NL)	Aneroid(AC)	Description
1820	1820NL	1820AC	Adult Complete Inflation System
1821	1821NL	1821AC	Child/Small Adult Complete Inflation System
1822	1822NL	1822AC	Infant Complete Inflation System
1823	1823NL	1823AC	Newborn Complete Inflation System
1824	1824NL	1824AC	Thigh Complete Inflation System
1825	1825NL	1825AC	Large Arm Complete Inflation System



1840

Inflation Bags

Cat. No.	Non Latex(NL)	Description
1840	1840NL	Adult Inflation Bag
1841	1841NL	Child/Small Adult Inflation Bag
1842	1842NL	Infant Inflation Bag
1844	1844NL	Thigh Inflation Bag
1845	1845NL	Large Arm Inflation Bag



1870

Calibrated® V-Lok® Cuffs

Cat. No.	Aneroid(AC)	Description
1870	1870AC	Adult Cuff Only
1871	1871AC	Child/Small Adult Cuff Only
1872	1872AC	Infant Cuff Only
1873		Newborn Cuff Only
1874	1874AC	Thigh Cuff Only
1868	1868AC	Large Arm Cuff Only



1890

Bulbs & Air-Flo® Control Valves

Cat. No.	Non Latex(NL)	Description
1890	1890NL	Bulb & Air-Flo® Control Valve
1891	1891NL	Bulb only
1892		Air-Flo® Control Valve only
1893	1893NL	Large Bulb only
1894	1894NL	Large Bulb & Air-Flo® Control Valve



2920



2910



2911

Latex Tubing and Connectors

Cat. No.	Non Latex(NL)	Description
2910	2910NL	Extendex® Tubing 8 ft. w/connectors
2919	2919NL	Extendex® Tubing 4 ft. w/connectors
2920		Tubing Connector, set
2911	2911NL	Straight Tubing, per foot

Note: Latex and Non-Latex inflation bags, bulbs and tubing are guaranteed for five years.

W.A. Baum Co. Inc.

SINCE 1916

ORIGINATOR AND MAKER OF ACCURATE BLOOD PRESSURE APPARATUS EXCLUSIVELY

Sound Investments

2700



The Baum Sprague Rappaport Type Stethoscope

Versatile - Equipped with modular adult and pediatric diaphragms, large, medium and small bell receivers which provide an extended range of auscultatory frequencies in one stethoscope.

Durable - Chrome plated zinc chestpiece and a rotatable chrome plated brass binaural.

Economical - Includes three types of eartips and a name tag/diaphragm replacement tool.

Cat. No.	Length	Non-Latex Tubing
2700	30 inches	Black
2701	30 inches	Gray
2702	30 inches	Blue

2720



The Baum Stainless Steel Dual Head Stethoscope

Multi-Functional - Designed for general health care applications. Built to our unique specifications.

Comfortable - Fitted with soft eartips and non-chill rings on the diaphragm and bell receivers.

Economical - Extra eartips and a replacement diaphragm are included.

Cat. No.	Length	Non-Latex Tubing
2720	31 inches	Black
2721	31 inches	Gray

2730



The Baum Stainless Steel Cardiology Stethoscope

High Performance - Delivers superior acoustic sensitivity. Built to our unique specifications.

Durable - Precision machined stainless steel dual head chestpiece, valve stem and binaural.

Comfortable - Fitted with soft eartips and non-chill rings on the diaphragm and bell receivers.

Cat. No.	Length	Non-Latex Tubing
2730	27 inches	Black
2731	27 inches	Gray

Guarantee: Baum Stethoscopes are guaranteed against defects in workmanship and material for a period of five years from date of purchase. During the guarantee period, stethoscope repair or replacement will be made without charge when returned, postage prepaid, unless it has been damaged in shipment, tampered with or abused.

W.A. Baum Co. Inc.

SINCE 1916

ORIGINATOR AND MAKER OF ACCURATE BLOOD PRESSURE APPARATUS EXCLUSIVELY