

Calibrated® V-Lok® Inflation System

Accurate: The Calibrated® V-Lok® System clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb.

Durable: Baum cuffs are made of tightly woven, heat set, Dacron™ polyester fabric that is double stitched at critical points for long wear. Hook and loop fasteners are evenly matched in holding strength to provide the optimal number of open-close cycles.

Comfortable: Cuffs conform to the shape of the limb and have no hard, stiff edges.

Latex-Free: Provide a high quality alternative for patients and practitioners who are hypersensitive to natural rubber latex. The inflation bulb and bag are a special blend of Neoprene that is resistant to attack from ozone, acids, bases, alcohols, and various hydrocarbons.

Safe: Each cuff is treated with an effective anti-microbial agent.

Cat. No.	Non Latex (NL)	Size	Calibrated® V-Lok®
1820	1820NL	Adult	Inflation System
1821	1821NL	Child/Sm. Adult	Inflation System
1822	1822NL	Infant	Inflation System
1823	1823NL	Newborn	Inflation System
1824	1824NL	Thigh	Inflation System
1825	1825NL	Large Arm	Inflation System
1870		Adult	Cuff Only
1871		Child/Sm. Adult	Cuff Only
1872		Infant	Cuff Only
1873		Newborn	Cuff Only
1874		Thigh	Cuff Only
1868		Large Arm	Cuff Only
1880	1880NL	Adult	Cuff & Inflation Bag
1881	1881NL	Child/Sm. Adult	Cuff & Inflation Bag
1882	1882NL	Infant	Cuff & Inflation Bag
1883	1883NL	Newborn	Cuff & Inflation Bag
1884	1884NL	Thigh	Cuff & Inflation Bag
1869	1869NL	Large Arm	Cuff & Inflation Bag



Aneroid Calibrated® V-Lok®

Accurate: The Calibrated® V-Lok® System clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb.

Durable: Baum cuffs are made of tightly woven, heat set, Dacron™ polyester fabric that is double stitched at critical points for long wear. Hook and loop fasteners are evenly matched in holding strength to provide the optimal number of open-close cycles.

Comfortable: Cuffs conform to the shape of the limb and have no hard, stiff edges.

Latex-Free: Provide a high quality alternative for patients and practitioners who are hypersensitive to natural rubber latex. The inflation bulb and bag are a special blend of Neoprene that is resistant to attack from ozone, acids, bases, alcohols, and various hydrocarbons.

Safe: Each cuff is treated with an effective anti-microbial agent.

Cat. No.	Non Latex (NL)	Size	Calibrated® V-Lok®
1820AC	1820ACNL	Adult	Inflation System
1821AC	1821ACNL	Child/Sm. Adult	Inflation System
1822AC	1822ACNL	Infant	Inflation System
1824AC	1824ACNL	Thigh	Inflation System
1825AC	1825ACNL	Large Arm	Inflation System
1870AC		Adult	Cuff Only
1871AC		Child/Sm. Adult	Cuff Only
1872AC		Infant	Cuff Only
1874AC		Thigh	Cuff Only
1868AC		Large Arm	Cuff Only
1880AC	1880ACNL	Adult	Cuff & Inflation Bag
1881AC	1881ACNL	Child/Sm. Adult	Cuff & Inflation Bag
1882AC	1882ACNL	Infant	Cuff & Inflation Bag
1884AC	1884ACNL	Thigh	Cuff & Inflation Bag
1869AC	1869ACNL	Large Arm	Cuff & Inflation Bag



Decorated & Disposable Cuffs

Accurate: The Calibrated® V-Lok® System clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb.

Decorated: Figures of animals and toys have been added in the hope that a child's anxiety and confusion may be diminished while the blood pressure is being measured.

Disposable Cuffs: are **Latex-Free** and fully interchangeable with our complete line of reusable cuffs and are intended for single patient use to help prevent cross contamination. They are available with two tubes (5220 or 5260 series) or one tube (5130 series). For use with manual sphygmomanometers or electronic non-invasive blood pressure monitors.

DECORATED CUFFS (D):

Cat. No.	Non Latex (NL)	Size	Calibrated® V-Lok®
1821D	1821DNL	Child/Sm. Adult	Inflation System
1822D	1822DNL	Infant	Inflation System
1881D	1881DNL	Child/Sm. Adult	Cuff & Inflation Bag
1882D	1882DNL	Infant	Cuff & Inflation Bag
1871D		Child/Sm. Adult	Cuff Only
1872D		Infant	Cuff Only

DOUBLE TUBE DISPOSABLE CUFFS:

Cat. No.	Size	Calibrated® V-Lok®	Packaging
5220	Adult	Inflation System	Box of 5
5221	Child/Sm. Adult	Inflation System	Box of 5
5222	Infant	Inflation System	Box of 5
5225	Large Arm	Inflation System	Box of 5
5260	Adult	Cuff & Inflation Bag	Box of 5
5261	Child/Sm. Adult	Cuff & Inflation Bag	Box of 5
5262	Infant	Cuff & Inflation Bag	Box of 5
5265	Large Arm	Cuff & Inflation Bag	Box of 5

SINGLE TUBE DISPOSABLE CUFFS:

Cat. No.	Size	Calibrated® V-Lok®	Packaging
5130	Adult	Cuff & Inflation Bag	Box of 5
5131	Child/Sm. Adult	Cuff & Inflation Bag	Box of 5
5132	Infant	Cuff & Inflation Bag	Box of 5
5135	Large Arm	Cuff & Inflation Bag	Box of 5

Available from your dealer.



Bulbs, Valves & Inflation Bags

Seamless: Baum bulbs and inflation bags are made of soft, smooth, dipped latex or Neoprene for seamless construction. They are protected against accelerated aging by means of effective anti-oxidants and anti-ozonants compounded in the liquid latex.

Latex-Free: Provide a high quality alternative for patients and practitioners who are hypersensitive to natural rubber latex. The Inflation bulb and bag are a special blend of Neoprene that is resistant to attack from ozone, acids, bases, alcohols, and various hydrocarbons.

Air-Flo® Control Valves

Instantly Responsive: Air entering the system is regulated by resilient valve tip of specially formulated silicone rubber designed for quick recovery and resistance to distortion and aging.

Smooth Operating: All key interior parts are stainless steel for corrosion-free, non-stick operation.

Durable: The Air-Flo® Control Valve body and air release mechanism are precision machined brass, nickel and chrome plated. Conical valve seat is integral with the valve body; it cannot loosen and cause leakage.

Cat. No.	Non Latex (NL)	Size	Calibrated® V-Lok®
Latex			
1840	1840NL	Adult	Inflation Bag
1841	1841NL	Child/Sm. Adult	Inflation Bag
1842	1842NL	Infant	Inflation Bag
1844	1844NL	Thigh	Inflation Bag
1845	1845NL	Large Arm	Inflation Bag

Cat. No.	Non Latex (NL)	Description
Latex		
1890	1890NL	Bulb & Air-Flo® Control Valve
1891	1891NL	Bulb only
1892		Air-Flo® Control Valve only
1893	1893NL	Large Bulb only
1894	1894NL	Large Bulb & Air-Flo® Control Valve

Baumanometer® Guarantee

Latex and Non-Latex inflation bags, bulbs and tubing are guaranteed for five years. This guarantee is invalid when other than genuine Baumanometer® parts are used.

Baumanometer, Standby, V-Lok, Extendex, Calibrated, Hi/Lo, Roll-By and Air-Flo are registered trademarks of W.A. Baum Co., Inc.



Calibrated® V-Lok® Inflation Systems and Components

W.A. Baum Co., Inc.
620 Oak Street
Copiague, NY 11726
www.wabaum.com

Toll Free: 1-888-281-6061
Tel: 631-226-3940
Fax: 631-226-3969
Email: info@wabaum.com

ISSC 8/05

PRINTED IN USA

W.A. Baum Co. Inc.

SINCE 1916

ORIGINATOR AND MAKER OF ACCURATE BLOOD PRESSURE APPARATUS EXCLUSIVELY

1880TA



Electronic NIBP Monitor Cuffs

Compatible: The two tube TA (Threaded Adapter) series fits connectors for Critikon, Kontnon and some Invivo monitors. The SSTA (Single Tube) series fits connectors for Welch-Allyn monitors. Luer Lock, Rectus and Colder Connectors (male and female) are available for a large variety of electronic NIBP monitors.

Durable: Baum cuffs are made of tightly woven, heat set, Dacron™ polyester fabric that is double stitched at critical points for long wear. Hook and loop fasteners are evenly matched in holding strength to provide the optimal number of open-close cycles.


Comfortable: Cuffs conform to the shape of the limb and have no hard, stiff edges.

Latex-Free: Provide a high quality alternative for patients and practitioners who are hypersensitive to natural rubber latex. The Inflation bulb and bag are a special blend of Neoprene that is resistant to attack from ozone, acids, bases, alcohols, and various hydrocarbons.

Safe: Each cuff is treated with an effective anti-microbial agent.

THREADED ADAPTERS (TA):

Cat. No.	Non Latex (NL)	Size	Calibrated® V-Lok®
1840TA	1840TANL	Adult	Inflation Bag
1841TA	1841TANL	Child/Sm. Adult	Inflation Bag
1842TA	1842TANL	Infant	Inflation Bag
1844TA	1844TANL	Thigh	Inflation Bag
1845TA	1845TANL	Large Arm	Inflation Bag
1880TA	1880TANL	Adult	Cuff & Inflation Bag
1881TA	1881TANL	Child/Sm. Adult	Cuff & Inflation Bag
1882TA	1882TANL	Infant	Cuff & Inflation Bag
1884TA	1884TANL	Thigh	Cuff & Inflation Bag
1869TA	1869TANL	Large Arm	Cuff & Inflation Bag
1840 SSTA	1840SSTANL	Adult	Inflation Bag
1841 SSTA	1841SSTANL	Child/Sm. Adult	Inflation Bag
1842 SSTA	1842SSTANL	Infant	Inflation Bag
1844 SSTA	1844SSTANL	Thigh	Inflation Bag
1845 SSTA	1845SSTANL	Large Arm	Inflation Bag
1880 SSTA	1880SSTANL	Adult	Cuff & Inflation Bag
1881 SSTA	1881SSTANL	Child/Sm. Adult	Cuff & Inflation Bag
1882 SSTA	1882SSTANL	Infant	Cuff & Inflation Bag
1884 SSTA	1884SSTANL	Thigh	Cuff & Inflation Bag
1869 SSTA	1869SSTANL	Large Arm	Cuff & Inflation Bag



Color Coded Cuff Sizes

All six cuff sizes are provided in Medium Blue unless otherwise specified. The color code for the cuffs shown above is as follows:

G - Gray **T** - Tan **DG** - Dark Green
MB - Medium Blue **S** - Scarlet **DB** - Dark Brown

Durable: Baum cuffs are made of tightly woven, heat set, Dacron™ polyester fabric that is double stitched at critical points for long wear. Hook and loop fasteners are evenly matched in holding strength to provide the optimal number of open-close cycles.

Comfortable: Cuffs conform to the shape of the limb and have no hard, stiff edges.

Safe: Cuffs are treated with an effective anti-microbial agent and can be paired with Non-Latex components for use with patients who are sensitive to natural rubber latex.

Quality: Baum components incorporate superior design, materials and workmanship to provide unsurpassed value.

Cat. No.	Non Latex (NL)	Size	Calibrated® V-Lok®
1820MB	1820MBNL	Adult	Inflation System
1821DG	1821DGNL	Child/Sm. Adult	Inflation System
1822T	1822TNL	Infant	Inflation System
1823G	1823GNL	Newborn	Inflation System
1824DB	1824DBNL	Thigh	Inflation System
1825S	1825SNL	Large Arm	Inflation System
1870MB		Adult	Cuff Only
1871DG		Child/Sm. Adult	Cuff Only
1872T		Infant	Cuff Only
1873G		Newborn	Cuff Only
1874DB		Thigh	Cuff Only
1868S		Large Arm	Cuff Only
1880MB	1880MBNL	Adult	Cuff & Inflation Bag
1881DG	1881DGNL	Child/Sm. Adult	Cuff & Inflation Bag
1882T	1882TNL	Infant	Cuff & Inflation Bag
1883G	1883GNL	Newborn	Cuff & Inflation Bag
1884DB	1884DBNL	Thigh	Cuff & Inflation Bag
1869S	1869SNL	Large Arm	Cuff & Inflation Bag

Available from your dealer.

Baumanometer® Guarantee

Latex and Non-Latex inflation bags, bulbs and tubing are guaranteed for five years. This guarantee is invalid when other than genuine Baumanometer® parts are used.

Baumanometer, Standby, V-Lok, Extendex, Calibrated, Hi/Lo, Roll-By and Air-Flo are registered trademarks of W.A. Baum Co., Inc.

W.A. Baum Co., Inc.

620 Oak Street

Copliague, NY 11726

www.wabaum.com

Toll Free: 1-888-281-6061

Tel: 631-226-3940

Fax: 631-226-3969

Email: info@wabaum.com