


One Standard. So Many Choices.



NIBP Cuff



Decorated Cuff



Latex-Free
(bulb and tubing)



Aneroid Cuff



- **Accurate:** The Calibrated® V-Lok® System clearly indicates whether or not the correct size cuff is being used.
- **Durable:** Baum cuffs are made of tightly woven, heat set, Dacron™ polyester fabric that is double stitched at critical points for long wear.
- **Comfortable:** Cuffs conform to the shape of the limb and have no hard, stiff edges. Each cuff is treated with an effective anti-microbial agent.
- **Latex-Free:** Provides a high quality alternative for those who are hypersensitive to natural rubber latex.
- **Guaranteed:** Latex and Non-Latex inflation bags, bulbs and tubing are guaranteed for five years.

For more information on the complete line of Baumanometer® Calibrated® V-Lok® Inflation Systems, call toll-free 1-888-281-6061 or email at info@wabaum.com



Bariatric Inflation System

W.A. Baum Co. Inc.

SINCE 1916

ORIGINATOR AND MAKER OF BLOOD PRESSURE APPARATUS EXCLUSIVELY