

COMPLETE CALIBRATED® V-LOK® INFLATION SYSTEMS

Standard	Qty	Cat. No.	Qty	Non-Latex	Size
		1820		1820NL	Adult
		1821		1821NL	Child/Sm. Adult
		1822		1822NL	Infant
		1823		1823NL	Newborn
		1824		1824NL	Thigh
		1825		1825NL	Large Arm




Aneroid (AC)	Qty	Cat. No.	Qty	Non-Latex	Size
		1820AC		1820ACNL	Adult
		1821AC		1821ACNL	Child/Sm. Adult
		1822AC		1822ACNL	Infant
		1824AC		1824ACNL	Thigh
		1825AC		1825ACNL	Large Arm



Decorated (D)	Qty	Cat. No.	Qty	Non-Latex	Size
		1821D		1821DNL	Child/Sm. Adult
		1822D		1822DNL	Infant



Bariatric	Qty	Cat. No.	Description
		1830	Bariatric Inflation System
		1830NL	Bariatric Non-Latex Inflation System

Range of Limb Circumferences  44-66 cm



CALIBRATED® V-LOK® CUFF AND INFLATION BAG SETS

Standard	Qty	Cat. No.	Qty	Non-Latex	Size
		1880		1880NL	Adult
		1881		1881NL	Child/Sm. Adult
		1882		1882NL	Infant
		1883		1883NL	Newborn
		1884		1884NL	Thigh
		1869		1869NL	Large Arm



Aneroid (AC)	Qty	Cat. No.	Qty	Non-Latex	Size
		1880AC		1880ACNL	Adult
		1881AC		1881ACNL	Child/Sm. Adult
		1882AC		1882ACNL	Infant
		1884AC		1884ACNL	Thigh
		1869AC		1869ACNL	Large Arm



Decorated (D)	Qty	Cat. No.	Qty	Non-Latex	Size
		1881D		1881DNL	Child/Sm. Adult
		1882D		1882DNL	Infant




CALIBRATED® V-LOK® CUFF AND INFLATION BAGS

CUFF SIZING CHART


Size	Fits Limb Circumference
Newborn	6 cm to 11 cm
Infant	10 cm to 19 cm
Child/Sm. Adult	18 cm to 26 cm
Adult	25 cm to 35 cm
Large Arm	33 cm to 47 cm
Thigh	46 cm to 66 cm

CALIBRATED® V-LOK® CUFFS

Standard	Qty	Cat. No.	Size
		1870	Adult
		1871	Child/Sm. Adult
		1872	Infant
		1873	Newborn
		1874	Thigh
		1868	Large Arm



Aneroid (AC)	Qty	Cat. No.	Size
		1870AC	Adult
		1871AC	Child/Sm. Adult
		1872AC	Infant
		1874AC	Thigh
		1868AC	Large Arm



Decorated (D)	Qty	Cat. No.	Size
		1871D	Child/Sm. Adult
		1872D	Infant



INFLATION BAGS

Double Tube	Qty	Cat. No.	Qty	Non-Latex	Size
		1840		1840NL	Adult
		1841		1841NL	Child/Sm. Adult
		1842		1842NL	Infant
		1843		1843NL	Newborn
		1844		1844NL	Thigh
		1845		1845NL	Large Arm



Single Tube	Qty	Cat. No.	Qty	Non-Latex	Size
		1840SS		1840SSNL	Adult
		1841SS		1841SSNL	Child/Sm. Adult
		1842SS		1842SSNL	Infant
		1843SS		1843SSNL	Newborn
		1844SS		1844SSNL	Thigh
		1845SS		1845SSNL	Large Arm





BAUM BULBS, AIR-FLO® CONTROL VALVES









Qty	Cat. No.	Qty	Non-Latex	Description
	1890		1890NL	Bulb & Valve
	1891		1891NL	Bulb only
	1892			Valve only
	1893		1893NL	Large Bulb only
	1894		1894NL	Large Bulb & Valve

LATEX & NON-LATEX TUBING


To order tubing for our various models, follow this guide: For Kompak, Desk, Hi/Lo® - 2 ft. straight; Standby® - 8 ft. Extendex®; Panel Wall, Wall Unit 33, Roll-By®, Wall Aneroid - 8 ft. Extendex®

Qty	Cat. No.	Qty	Non-Latex	Description
	2910		2910NL	Extendex® Tubing (8 ft)
	2911		2911NL	Straight Tubing, per ft.
	2900			Straight Tubing, 100 ft. Reel
	2950		2950NL	Straight Tubing, 50 ft. Reel



TUBING CONNECTORS

Qty	Cat. No.	Description
	2920	Tubing Connectors, set (#2921 and #2922)
	2921	Tubing Connector, male only
	2922	Tubing Connector, female only
	2923	Tubing Connector, for Extendex®, male only
	2934	Straight Connector
	2935	Connector, Extendex®-to-Straight Tubing



BAUMANOMETER® MYLAR® CLAD GLASS CARTRIDGE TUBES

Qty	Cat. No.	Description
	2100	260 mmHg for Kompak, Hi/Lo®
	2200	300 mmHg for Standby®
	2300	300 mmHg for Desk, Panel Wall, WU 33 Models

BAUM ACCESSORIES

Qty	Cat. No.	Description
	2420	Nylon Coated Steel Wall Basket
	2422	Large Nylon Coated Steel Wall Basket
	2421	Stainless Steel Wall Clip

TEST KITS

Qty	Cat. No.	Qty	Non-Latex	Description
	2940		2940NL	Calibration "Y" Connector
	2941		2941NL	Above with Bulb

NIBP MONITOR CUFFS & INFLATION BAGS

Cuff & Bag w/2 Threaded Adapters (Double Tube)

Qty	Cat. No.	Qty	Non-Latex	Size
	1880TA		1880TANL	Adult
	1881TA		1881TANL	Child/Small Adult
	1882TA		1882TANL	Infant
	1884TA		1884TANL	Thigh
	1869TA		1869TANL	Large Arm

Cuff & Bag w/1 Threaded Adapter (Single Tube)

Qty	Cat. No.	Qty	Non-Latex	Size
	1880SSTA		1880SSTANL	Adult
	1881SSTA		1881SSTANL	Child/Small Adult
	1882SSTA		1882SSTANL	Infant
	1884SSTA		1884SSTANL	Thigh
	1869SSTA		1869SSTANL	Large Arm

Cuff & Bag w/1 Female Luer Lock (Single Tube)

Qty	Cat. No.	Qty	Non-Latex	Size
	1881SSLL		1881SSLLNL	Child/Small Adult
	1882SSLL		1882SSLLNL	Infant
	1884SSLL		1884SSLLNL	Thigh
	1869SSLL		1869SSLLNL	Large Arm

Cuff & Inflation Bag w/1/Male Rectus Connector (Single Tube)

Qty	Cat. No.	Qty	Non-Latex	Size
	1880SSOC		1880SSOCNL	Adult
	1881SSOC		1881SSOCNL	Child/Small Adult
	1882SSOC		1882SSOCNL	Infant
	1884SSOC		1884SSOCNL	Thigh
	1869SSOC		1869SSOCNL	Large Arm

Cuff & Inflation Bag w/ 1 Male & 1 Female Colder Connector (Double Tube)

Qty	Cat. No.	Qty	Non-Latex	Size
	1880MFCC		1880MFCCNL	Adult
	1881MFCC		1881MFCCNL	Child/Small Adult
	1882MFCC		1882MFCCNL	Infant
	1884MFCC		1884MFCCNL	Thigh
	1869MFCC		1869MFCCNL	Large Arm

Cuff & Inflation Bag w/ 2 Female Colder Connectors (Double Tube)

Qty	Cat. No.	Qty	Non-Latex	Size
	1880FFCC		1880FFCCNL	Adult
	1881FFCC		1881FFCCNL	Child/Small Adult
	1882FFCC		1882FFCCNL	Infant
	1884FFCC		1884FFCCNL	Thigh
	1869FFCC		1869FFCCNL	Large Arm

Cuff & Inflation Bag w/ 1 Female Colder Connector (Single Tube)


Qty	Cat. No.	Qty	Non-Latex	Size
	1880SSCC		1880SSCCNL	Adult
	1881SSCC		1881SSCCNL	Child/Small Adult
	1882SSCC		1882SSCCNL	Infant
	1884SSCC		1884SSCCNL	Thigh
	1869SSCC		1869SSCCNL	Large Arm

Available from your dealer.

Baumanometer® Guarantee

Latex and Non-Latex inflation bags, bulbs and tubing are guaranteed for five years. This guarantee is invalid when other than genuine Baumanometer® parts are used.

Baumanometer, Standby, V-Lok, Extendex, Calibrated, Hi-Lo, Roll-By and Air-Flo are registered trademarks of W.A. Baum Co., Inc.



Baum Replacement Parts

W.A. Baum Co., Inc.
620 Oak Street
Copiague, NY 11726
www.wabaum.com

Toll Free: 1-888-281-6061
Tel: 631-226-3940
Fax: 631-226-3969
Email: info@wabaum.com

W.A. Baum Co. Inc.

SINCE 1916
ORIGINATOR AND MAKER OF ACCURATE
BLOOD PRESSURE APPARATUS EXCLUSIVELY