

Available from your distributor:



### DESK ANEROID

- 300 mmHg Aneroid instrument mounted at a 20° angle for ease of viewing.
- Complies with ANSI/AAMI SP-10 for accuracy and performance (1% full scale  $\pm$  3 mmHg).
- Steel stand finished with Neutral Gray textured baked enamel and non-skid base..
- Neutral White ABS plastic manometer with chrome plated steel bracket.
- Unique high contrast scale and red pointer increase visibility and reduce parallax error.
- Manufactured with a low-friction precision movement and twin pre-seasoned phosphor-bronze diaphragm capsules.
- Available with Latex or **Non-Latex** components.
- **Calibrated® V-Lok® Inflation System.**  
**Accurate:** The Calibrated® V-Lok® system clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb. **Comfortable:** Cuffs conform to the shape of the limb and have no hard, stiff edges. **Durable:** Baum cuffs are made of tightly woven, heat set Dacron™ polyester fabric. They are double stitched at critical points for long wear and treated with an effective anti-microbial agent.  
**Non-Latex:** Available for those who are hypersensitive to natural rubber latex.

**Instrument Finish:** Neutral White ABS plastic  
**Steel Stand Finish:** Neutral Gray baked enamel with non-skid base.  
**Standard Cuff Color:** Medium Blue

#### DESK ANEROID with 25" Straight Tubing and Steel Stand:

Cat. No.	Non-Latex (NL)	Cuff Size	Inflation System
0920	0920NL	Adult	Calibrated® V-Lok®
0921	0921NL	Child/Sm. Adult	Calibrated® V-Lok®
0922	0922NL	Infant	Calibrated® V-Lok®
0923	0923NL	Newborn	Calibrated® V-Lok®
0924	0924NL	Thigh	Calibrated® V-Lok®
0925	0925NL	Large Arm	Calibrated® V-Lok®



### Clinical Sphygmomanometers

### Baum Aneroid Guarantee

Every Baum Aneroid instrument is guaranteed against defects in workmanship and material for a period of one year from date of purchase. This guarantee is invalid when other than genuine W. A. Baum parts are used.

We recommend that every aneroid instrument be checked for accuracy periodically by following our "Testing Instrument Accuracy" procedure ( check our website ) to determine if it needs adjustment.

Adjustments will be made to Baum Aneroid instruments without charge when returned postage prepaid, unless the instrument has been damaged in shipment, tampered with or abused.

**W. A. Baum Co., Inc.**  
 620 Oak Street  
 Copiague, NY 11726  
[www.wabaum.com](http://www.wabaum.com)

Toll Free: 1-888-281-6061  
 Tel: 631-226-3940  
 Fax: 631-226-3969  
 Email: [info@wabaum.com](mailto:info@wabaum.com)

**W.A. Baum Co. Inc.**

SINCE 1916  
 ORIGINATOR AND MAKER OF ACCURATE  
 BLOOD PRESSURE APPARATUS EXCLUSIVELY

**0950****WALL ANEROID**

- 300 mmHg swivel-type (110° arc) clinical sphygmomanometer.
- Complies with ANSI/AAMI SP-10 for accuracy and performance (1% full scale  $\pm$  3 mmHg).
- Neutral White ABS plastic manometer with chrome plated steel bracket.
- Unique high contrast scale and red pointer.
- Available with Latex or **Non-Latex** components.
- Available with large or small basket.
- **Calibrated® V-Lok® Inflation System.**

**Accurate:** The Calibrated® V-Lok® system clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb. **Comfortable:** Cuffs conform to the shape of the limb and have no hard, stiff edges. **Durable:** Baum cuffs are made of tightly woven, heat set Dacron™ polyester fabric. They are double stitched at critical points for long wear and treated with an effective anti-microbial agent.

**Non-Latex:** Available for those who are hypersensitive to natural rubber latex.

**Instrument Finish:** Neutral White ABS Plastic

**Standard Cuff Color:** Medium Blue

**Basket Color:** Neutral White (standard)

**WALL ANEROID** with 8' Extendex® Tubing and Wall Basket:

Cat. No.	Non-Latex (NL)	Cuff Size	Inflation System
0950	0950NL	Adult	Calibrated® V-Lok®
0951	0951NL	Child/Sm. Adult	Calibrated® V-Lok®
0952	0952NL	Infant	Calibrated® V-Lok®
0953	0953NL	Newborn	Calibrated® V-Lok®
0954	0954NL	Thigh	Calibrated® V-Lok®
0955	0955NL	Large Arm	Calibrated® V-Lok®
0960	0960NL	Adult (No Basket)	Calibrated® V-Lok®
0933	0933NL	None (w/8 ft. Extendex® tubing only)	No Basket.
0970	0970NL	Adult, Child/Sm. Adult & Large Arm w/8 ft. Extendex® tubing, large storage basket (#2422)	

**Wall Instrument Accessories**

Cat. No.	Description
2420	Wall Basket (Nylon Coated Steel)
2422	Large Wall Basket (Nylon Coated Steel)
2421	Wall Clip (Stainless Steel)

**1050****POCKET ANEROID**

- 300 mmHg clinical sphygmomanometer.
- Complies with ANSI/AAMI SP-10 for accuracy and performance (1% full scale  $\pm$  3 mmHg).
- Diecast zinc manometer housing.
- Precision low friction movement with twin pre-seasoned phosphor-bronze diaphragm capsules.
- Unique high contrast scale and red pointer.
- Available with Latex or **Non-Latex** components.
- **Calibrated® V-Lok® Inflation System.**

**Accurate:** The Calibrated® V-Lok® system clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb. **Comfortable:** Cuffs conform to the shape of the limb and have no hard, stiff edges. **Durable:** Baum cuffs are made of tightly woven, heat set Dacron™ polyester fabric. They are double stitched at critical points for long wear and treated with an effective anti-microbial agent.

**Non-Latex:** Available for those who are hypersensitive to natural rubber latex.

**Finish:** Warm Gray textured baked enamel

**Standard Cuff Color:** Medium Blue

**POCKET ANEROID** with Calibrated® V-Lok® and Zippered Case:

Cat. No.	Non-Latex (NL)	Cuff Size	Inflation System
1050	1050NL	Adult	Calibrated® V-Lok®
1051	1051NL	Child/Sm. Adult	Calibrated® V-Lok®
1052	1052NL	Infant	Calibrated® V-Lok®
1054	1054NL	Thigh (No Case)	Calibrated® V-Lok®
1055	1055NL	Large Arm	Calibrated® V-Lok®

**Baum Aneroid Guarantee**

Every Baum Aneroid instrument is guaranteed against defects in workmanship and material for a period of one year from date of purchase. This guarantee is invalid when other than genuine W. A. Baum parts are used.

We recommend that every aneroid instrument be checked for accuracy periodically by following our "Testing Instrument Accuracy" procedure (check our website) to determine if it needs adjustment.

Adjustments will be made to Baum Aneroid instruments without charge when returned postage prepaid, unless the instrument has been damaged in shipment, tampered with or abused.

**1150****ROLL-BY®**

- 300 mmHg clinical sphygmomanometer.
- Complies with ANSI/AAMI SP-10 for accuracy and performance (1% full scale  $\pm$  3 mmHg).
- Unique, low reflection, high contrast scale and red pointer.
- Diecast aluminum 5-wheel base (20" dia.) with 2" hooded casters and 2 lb. cast iron weight for increased stability.
- Manometer housing is mounted at a 20° angle for ease of viewing.
- Heavy gauge machined aluminum tubing.
- Available with Latex or **Non-Latex** components.
- **Calibrated® V-Lok® Inflation System.**

**Accurate:** The Calibrated® V-Lok® system clearly indicates whether or not the correct size cuff is being used. The center of the bladder is clearly marked on the inner and outer surfaces of the cuff for correct application on either the right or left limb. **Comfortable:** Cuffs conform to the shape of the limb and have no hard, stiff edges. **Durable:** Baum cuffs are made of tightly woven, heat set Dacron™ polyester fabric. They are double stitched at critical points for long wear and treated with an effective anti-microbial agent.

**Non-Latex:** Available for those who are hypersensitive to natural rubber latex.

**Finish:** Neutral White ABS plastic

**Standard Cuff Color:** Medium Blue

**ROLL-BY®** with 8' Extendex® Tubing and Large Circular Basket:

Cat. No.	Non-Latex (NL)	Cuff Size(s)	Inflation System
1150	1150NL	Adult	Calibrated® V-Lok®
1151	1151NL	Child/Sm. Adult	Calibrated® V-Lok®
1152	1152NL	Infant	Calibrated® V-Lok®
1153	1153NL	Newborn	Calibrated® V-Lok®
1154	1154NL	Thigh	Calibrated® V-Lok®
1155	1155NL	Large Arm	Calibrated® V-Lok®
1170	1170NL	Adult, Child, and Large Arm	Calibrated® V-Lok®
1140	1140NL	Roll-By® Manometer Only (Models 1150/1170)	